



DIRECT HITS

KEYS TO SUCCESSFUL SELLING

Standard Industrial Classification
 3841– Surgical and medical instruments and apparatus

Application
 Active cooling devices have separate mechanisms (e.g., pumps) that attach to garments, circulating coolant through tubes that require quick and clean connections.

In This Application
 BreakAway™ couplings are ideal for medical cooling device connections because they are designed to disconnect when pull force is applied without requiring the user to depress a latching mechanism. This allows a quick, leak-free break that prevents messy spills and equipment damage.

Media
 Water and coolants

Pressure
 0 - 20 psi

Temperature
 0° - 150° F

- Colder Part Numbers**
- BACD16004 - Acetal panel mount body with shutoff and Buna-N o-ring
 - BACD22004 - Acetal in-line insert with shutoff and Buna-N o-ring
- Note: Available in black or white acetal

- Features & Benefits**
- 1. Reliable break away design**
 - Rapid automatic disconnection
 - Easy hands-free operation
 - Dependable for multiple disconnections
 - 2. Automatic dual shutoff valves**
 - Avoid leaks and spills
 - Reduce clean-up time
 - Safe and clean fluid handling

Active Cooling Device



Surgeons and patients use liquid-cooled vests to keep comfortable during treatment and/or long operations.



- Features & Benefits**
- 3. Streamlined, flush finish**
 - Easy-to-use for maximum productivity
 - Intuitive compact component
 - Developed from proven platform



All sales are subject to Colder Products Company's limited express warranty set forth in the Colder catalog. Contact your local distributor or Colder Customer Service for warranty provisions.
 1001 Westgate Drive, St. Paul, MN 55114 U.S.A. • Phone 651-645-0091 • Fax 651-645-5404 • www.colder.com

SIC
3841

Direct Hit
BAC1